

Product Datasheet

Human N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit (orb1205937)

Description	Human N-Acylsphingosine Amidohydrolase 1 (ASAH1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13510
Sensitivity	0.61ng/mL
Entrez	427
Expiration Date	Please enquire.