

Product Datasheet

Human Actin Related Protein 2/3 Complex Subunit 2 (ARPC2) ELISA Kit (orb1205944)

Description	Human Actin Related Protein 2/3 Complex Subunit 2 (ARPC2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	O15144
Sensitivity	0.068ng/mL
Entrez	10109
Expiration Date	Please enquire.