

Product Datasheet

Human Amyloid Beta Precursor Like Protein 1 (APLP1) ELISA Kit (orb1205961)

Description	Human Amyloid Beta Precursor Like Protein 1 (APLP1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P51693
Sensitivity	0.118ng/mL
Entrez	333
Expiration Date	Please enquire.