

Product Datasheet

Human Amine Oxidase, Copper Containing 3 (AOC3) ELISA Kit (orb1205968)

Description	Human Amine Oxidase, Copper Containing 3 (AOC3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q16853
Sensitivity	0.051ng/mL
Entrez	8639
Expiration Date	Please enquire.