

## Product Datasheet

### Human Acyloxyacyl Hydrolase (AOAH) ELISA Kit (orb1205969)

|                          |                                                                                |
|--------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>       | Human Acyloxyacyl Hydrolase<br>(AOAH) ELISA Kit                                |
| <b>Reactivity</b>        | Human                                                                          |
| <b>Range</b>             | 0.156-10ng/mL                                                                  |
| <b>Note</b>              | For research use only                                                          |
| <b>Application notes</b> | sandwich                                                                       |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates,<br>cell lysates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>P28039</b>                                                                  |
| <b>Sensitivity</b>       | 0.059ng/mL                                                                     |
| <b>Entrez</b>            | <b>313</b>                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                                |