

Product Datasheet

Human Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) ELISA Kit (orb1205985)

Description	Human Alkaline Phosphatase, Placental Like Protein 2 (ALPPL2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P10696
Sensitivity	0.113ng/mL
Entrez	251
Expiration Date	Please enquire.