

Product Datasheet

Human Arachidonate Lipoxygenase 3 (ALOXE3) ELISA Kit (orb1205986)

Description	Human Arachidonate Lipoxygenase 3 (ALOXE3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9BYJ1
Sensitivity	0.114ng/mL
Entrez	59344
Expiration Date	Please enquire.