

Product Datasheet

Human Aldehyde Dehydrogenase, Mitochondrial (ALDM) ELISA Kit (orb1205989)

Description	Human Aldehyde Dehydrogenase, Mitochondrial (ALDM) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P05091
Sensitivity	0.078ng/mL
Entrez	217
Expiration Date	Please enquire.