

Product Datasheet

Human Aminolevulinate Delta Dehydratase (ALAD) ELISA Kit (orb1205995)

Description	Human Aminolevulinate Delta Dehydratase (ALAD) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, erythrocyte lysates or other biological fluids.
Uniprot ID	P13716
Sensitivity	0.128ng/mL
Entrez	210
Expiration Date	Please enquire.