

Product Datasheet

Human Aryl Hydrocarbon Receptor Repressor (AHRR) ELISA Kit (orb1206000)

Description	Human Aryl Hydrocarbon Receptor Repressor (AHRR) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	A9YTQ3
Sensitivity	0.062ng/mL
Entrez	57491
Expiration Date	Please enquire.