

Product Datasheet

Human Anterior Gradient Protein 2 (AGR2) ELISA Kit (orb1206003)

Description	Human Anterior Gradient Protein 2 (AGR2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O95994
Sensitivity	0.061ng/mL
Entrez	10551
Expiration Date	Please enquire.