

Product Datasheet

Human Adrenergic Receptor Beta Kinase 1 (ADRBK1) ELISA Kit (orb1206009)

Description	Human Adrenergic Receptor Beta Kinase 1 (ADRBK1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P25098
Sensitivity	0.105ng/mL
Entrez	156
Expiration Date	Please enquire.