

Product Datasheet

Human Adrenergic Receptor Beta 1 (ADRB1) ELISA Kit (orb1206012)

Description	Human Adrenergic Receptor Beta 1 (ADRB1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P08588
Sensitivity	0.058ng/mL
Entrez	153
Expiration Date	Please enquire.