

Product Datasheet

Human Adrenergic Receptor Alpha 2C (ADRa2C) ELISA Kit (orb1206013)

Description	Human Adrenergic Receptor Alpha 2C (ADRa2C) ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P18825
Sensitivity	0.237ng/mL
Entrez	152
Expiration Date	Please enquire.