

Product Datasheet

Human A Disintegrin And Metalloprotease 28 (ADAM28) ELISA Kit (orb1206025)

Description	Human A Disintegrin And Metalloprotease 28 (ADAM28) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9UKQ2
Sensitivity	0.061ng/mL
Entrez	10863
Expiration Date	Please enquire.