

Product Datasheet

Human A Disintegrin And Metalloprotease 19 (ADAM19) ELISA Kit (orb1206027)

Description	Human A Disintegrin And Metalloprotease 19 (ADAM19) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9H013
Sensitivity	0.108ng/mL
Entrez	8728
Expiration Date	Please enquire.