

Product Datasheet

Human ATP Binding Cassette Transporter C6 (ABCC6) ELISA Kit (orb1206040)

Description	Human ATP Binding Cassette Transporter C6 (ABCC6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O95255
Sensitivity	0.126ng/mL
Entrez	368
Expiration Date	Please enquire.