

Product Datasheet

Human ATP Binding Cassette Transporter C4 (ABCC4) ELISA Kit (orb1206042)

Description	Human ATP Binding Cassette Transporter C4 (ABCC4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O15439
Sensitivity	0.118ng/mL
Entrez	10257
Expiration Date	Please enquire.