

Product Datasheet

Human ATP Binding Cassette Transporter C13 (ABCC13) ELISA Kit (orb1206045)

Description	Human ATP Binding Cassette Transporter C13 (ABCC13) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9NSE7
Sensitivity	0.057ng/mL
Entrez	150000
Expiration Date	Please enquire.