

Product Datasheet

Human ATP Binding Cassette Transporter C12 (ABCC12) ELISA Kit (orb1206046)

Description	Human ATP Binding Cassette Transporter C12 (ABCC12) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q96J65
Sensitivity	0.072ng/mL
Entrez	94160
Expiration Date	Please enquire.