

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily, Member 4 (TNFRSF4) ELISA Kit (orb1206054)

Description	Human Tumor Necrosis Factor Receptor Superfamily, Member 4 (TNFRSF4) ELISA Kit
Reactivity	Human
Range	7.81-500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P43489
Sensitivity	2.89pg/mL
Entrez	7293
Expiration Date	Please enquire.