

## Product Datasheet

### Bovine Pleiotrophin (PTN) ELISA Kit (orb1206064)

|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Description</b>       | Bovine Pleiotrophin (PTN)<br>ELISA Kit       |
| <b>Reactivity</b>        | Bovine                                       |
| <b>Range</b>             | 31.2-2000pg/mL                               |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | sandwich                                     |
| <b>Sample Types</b>      | serum, plasma or other<br>biological fluids. |
| <b>Uniprot ID</b>        | <b>P21782</b>                                |
| <b>Sensitivity</b>       | 11.4pg/mL                                    |
| <b>Entrez</b>            | <b>280904</b>                                |
| <b>Expiration Date</b>   | Please enquire.                              |