

Product Datasheet

Mouse Tumor Necrosis Factor Alpha Induced Protein 6 (TNFaIP6) ELISA Kit (orb1206065)

Description	Mouse Tumor Necrosis Factor Alpha Induced Protein 6 (TNFaIP6) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	008859
Sensitivity	0.054ng/mL
Entrez	21930
Expiration Date	Please enquire.