

Product Datasheet

Rabbit Vascular Cell Adhesion Molecule 1 (VCAM1) ELISA Kit (orb1206095)

Description	Rabbit Vascular Cell Adhesion Molecule 1 (VCAM1) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	G1T1S9
Sensitivity	5.4pg/mL
Entrez	100008901
Expiration Date	Please enquire.