

Product Datasheet

Rat Visceral Adipose Tissue Derived Serine Protease Inhibitor (Vaspin) ELISA Kit (orb1206097)

Description	Rat Visceral Adipose Tissue Derived Serine Protease Inhibitor (Vaspin) ELISA Kit
Reactivity	Rat
Range	93.7-6000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q8R4Z1
Sensitivity	33.7pg/mL
Entrez	191570
Expiration Date	Please enquire.