

Product Datasheet

Rabbit Tissue Inhibitors Of Metalloproteinase 1 (TIMP1) ELISA Kit (orb1206132)

Description	Rabbit Tissue Inhibitors Of Metalloproteinase 1 (TIMP1) ELISA Kit
Reactivity	Rabbit
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P20614
Sensitivity	0.119ng/mL
Entrez	100009047
Expiration Date	Please enquire.