

Product Datasheet

Mouse Signal Transducer And Activator Of Transcription 3 (STAT3) ELISA Kit (orb1206152)

Description	Mouse Signal Transducer And Activator Of Transcription 3 (STAT3) ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42227
Sensitivity	0.31ng/mL
Entrez	20848
Expiration Date	Please enquire.