

Product Datasheet

Rat Sex Determining Region Y Box Protein 2 (SOX2) ELISA Kit (orb1206159)

Description	Rat Sex Determining Region Y Box Protein 2 (SOX2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	D4A543
Sensitivity	0.064ng/mL
Entrez	499593
Expiration Date	Please enquire.