

Product Datasheet

Rat Superoxide Dismutase, Copper Chaperone (SOD4) ELISA Kit (orb1206160)

Description	Rat Superoxide Dismutase, Copper Chaperone (SOD4) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9JK72
Sensitivity	0.064ng/mL
Entrez	84485
Expiration Date	Please enquire.