

Product Datasheet

Rat Mothers Against Decapentaplegic Homolog 7 (Smad7) ELISA Kit (orb1206166)

Description	Rat Mothers Against Decapentaplegic Homolog 7 (Smad7) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O88406
Sensitivity	0.055ng/mL
Entrez	81516
Expiration Date	Please enquire.