

Product Datasheet

Rat Selectin, Leukocyte (SELL) ELISA Kit (orb1206174)

Description	Rat Selectin, Leukocyte (SELL) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P30836
Sensitivity	0.041ng/mL
Entrez	29259
Expiration Date	Please enquire.