

Product Datasheet

Rabbit S100 Calcium Binding Protein B (S100B) ELISA Kit (orb1206182)

Description	Rabbit S100 Calcium Binding Protein B (S100B) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6YNR6
Sensitivity	4.5pg/mL
Entrez	100009495
Expiration Date	Please enquire.