

Product Datasheet

Rat Peroxisome Proliferator Activated Receptor Gamma (PPAR γ) ELISA Kit (orb1206221)

Description	Rat Peroxisome Proliferator Activated Receptor Gamma (PPAR γ) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O88275
Sensitivity	0.058ng/mL
Entrez	25664
Expiration Date	Please enquire.