

Product Datasheet

Rat Protease Activated Receptor 4 (PAR4) ELISA Kit (orb1206258)

Description	Rat Protease Activated Receptor 4 (PAR4) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q920E0
Sensitivity	0.059ng/mL
Entrez	116498
Expiration Date	Please enquire.