

## Product Datasheet

# Rat N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit (orb1206276)

|                          |                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Rat N-Terminal Pro Brain<br>Natriuretic Peptide (NT-ProBNP)<br>ELISA Kit                                      |
| <b>Reactivity</b>        | Rat                                                                                                           |
| <b>Range</b>             | 78.1-5000pg/mL                                                                                                |
| <b>Note</b>              | For research use only                                                                                         |
| <b>Application notes</b> | sandwich                                                                                                      |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates, cell lysates, cell<br>culture supernates or other<br>biological fluids. |
| <b>Sensitivity</b>       | 27pg/mL                                                                                                       |
| <b>Expiration Date</b>   | Please enquire.                                                                                               |