

Product Datasheet

Mouse Mannose Receptor C Type 1 (MRC1) ELISA Kit (orb1206319)

Description	Mouse Mannose Receptor C Type 1 (MRC1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q61830
Sensitivity	0.059ng/mL
Entrez	17533
Expiration Date	Please enquire.