

Product Datasheet

Rabbit Matrix Metalloproteinase 3 (MMP3) ELISA Kit (orb1206327)

Description	Rabbit Matrix Metalloproteinase 3 (MMP3) ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P28863
Sensitivity	0.061ng/mL
Entrez	100009111
Expiration Date	Please enquire.