

Product Datasheet

Rabbit Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit (orb1206343)

Description	Rabbit Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P46632
Sensitivity	6.2pg/mL
Entrez	100008984
Expiration Date	Please enquire.