

Product Datasheet

Mouse Mucosal Addressin Cell Adhesion Molecule 1 (MAdCAM1) ELISA Kit (orb1206364)

Description	Mouse Mucosal Addressin Cell Adhesion Molecule 1 (MAdCAM1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, urine or other biological fluids.
Uniprot ID	Q61826
Sensitivity	0.059ng/mL
Entrez	17123
Expiration Date	Please enquire.