

Product Datasheet

Mouse Lectin Galactoside Binding, Soluble 3 Binding Protein (LGALS3BP) ELISA Kit (orb1206374)

Description	Mouse Lectin Galactoside Binding, Soluble 3 Binding Protein (LGALS3BP) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q07797
Sensitivity	0.125ng/mL
Entrez	19039
Expiration Date	Please enquire.