

Product Datasheet

Mouse Keratin 18 (KRT18) ELISA Kit (orb1206387)

Description	Mouse Keratin 18 (KRT18) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P05784
Sensitivity	25pg/mL
Entrez	16668
Expiration Date	Please enquire.