

Product Datasheet

Rat Heat Shock 60kD Protein 1, Chaperonin (HSPD1) ELISA Kit (orb1206419)

Description	Rat Heat Shock 60kD Protein 1, Chaperonin (HSPD1) ELISA Kit
Reactivity	Rat
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P63039
Sensitivity	1.27ng/mL
Entrez	63868
Expiration Date	Please enquire.