

Product Datasheet

Mouse Heat Shock 10kDa Protein 1 (HSP10) ELISA Kit (orb1206425)

Description	Mouse Heat Shock 10kDa Protein 1 (HSP10) ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q64433
Sensitivity	0.65ng/mL
Entrez	15528
Expiration Date	Please enquire.