

Product Datasheet

Rat Hexosaminidase B Beta (HEXb) ELISA Kit (orb1206435)

Description	Rat Hexosaminidase B Beta (HEXb) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q6AXR4
Sensitivity	0.59ng/mL
Entrez	294673
Expiration Date	Please enquire.