

Product Datasheet

Rat Glutathione S Transferase Alpha 3 (GSTa3) ELISA Kit (orb1206446)

Description	Rat Glutathione S Transferase Alpha 3 (GSTa3) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P46418
Sensitivity	0.079ng/mL
Entrez	494500
Expiration Date	Please enquire.