

Product Datasheet

Rat Glycogen Synthase Kinase 3 Alpha (GSK3a) ELISA Kit (orb1206448)

Description	Rat Glycogen Synthase Kinase 3 Alpha (GSK3a) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P18265
Sensitivity	0.227ng/mL
Entrez	50686
Expiration Date	Please enquire.