

Product Datasheet

Rat Gonadotropin Releasing Hormone (GnRH) ELISA Kit (orb1206458)

Description	Rat Gonadotropin Releasing Hormone (GnRH) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07490
Sensitivity	29pg/mL
Entrez	25194
Expiration Date	Please enquire.