

Product Datasheet

Mouse Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH) ELISA Kit (orb1206474)

Description	Mouse Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P16858
Sensitivity	0.055ng/mL
Entrez	14433
Expiration Date	Please enquire.