

Product Datasheet

Mouse Receptor II For The Fc Region Of Immunoglobulin E (FceRII) ELISA Kit (orb1206493)

Description	Mouse Receptor II For The Fc Region Of Immunoglobulin E (FceRII) ELISA Kit
Reactivity	Mouse
Range	1.56-100pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P20693
Sensitivity	0.66pg/mL
Entrez	14128
Expiration Date	Please enquire.