

Product Datasheet

Mouse Fatty Acid Binding Protein 5, Epidermal (FABP5) ELISA Kit (orb1206497)

Description	Mouse Fatty Acid Binding Protein 5, Epidermal (FABP5) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, breast milk or other biological fluids.
Uniprot ID	Q05816
Sensitivity	0.054ng/mL
Entrez	16592
Expiration Date	Please enquire.